

Goal: Develop a system that enables 19 or fewer vials to be utilized for protocol development.

Theory: A small batch will freeze dry faster than a large batch in a lab freeze dryer.

First Experimental Set: A full tray in a lab freeze dryer and then run 19 vials in the same dryer using the same cycle (using CM-Pirani convergence to indicate end of primary drying).

Experimental Group # 1 Water

Full tray in REVO Vs 19 Vials in REVO

Run a REVO with a full tray of vials and time
 Experimental Design
 3ml fill in a 10 ml vial
 Water
 Freezing step ramp 0.5C/min to -40C, hold 60 mins
 Primary drying shelf set point +20C
 Vacuum set point 100 mT

Run a REVO with a 19 vials and time
 Experimental Design
 3ml fill in a 10 ml vial
 Water
 Freezing step ramp 0.5C/min to -40C, hold 60 mins
 Primary drying shelf set point +20C
 Vacuum set point 100 mT


REVO Full Tray
636 minutes in primary drying

636 minutes is the target for all other water runs

REVO 19 Vials
512 minutes in primary drying

Results

19 vials dry quicker than a full tray based on end of primary drying being detected by CM vs Pirani convergence



Results: 19 Vials Take less time. Why?
Assumption: Radiant energy has a significant effect on the freeze drying process via edge vial effect. Smaller batch runs are more susceptible to this effect because a larger percentage of the vials are "edge effect vials".

Experimental Group # 1

5% Sucrose

Full tray in REVO

Run a REVO with a full tray of vials and time

Experimental Design
2ml fill in a 10 ml vial
5% sucrose by weight
Freezing step -10 hold 30 mins, -40C hold 120 mins
Primary drying shelf set point -20C
Vacuum set point 60 mT




REVO Full Tray
2026 minutes in primary drying

2026 minutes is the target for all other sucrose runs

Results

Full tray took 2026 mins based on end of primary drying being detected by CM vs Pirani convergence



Experimental Group #2

Water


Eliminate radiant energy to minimize edge effect vials by cooling and/or shielding the chamber walls of the freeze dryer.

Results

The minimization of radiant energy does not eliminate edge effect.

Theory: Eliminate radiant energy, and thus edge effect vials, and your small batch runs will be more like the larger batch runs.

Second Experimental Set: Eliminate radiant energy from the freeze dryer when using 19 vials. Build a small freeze dryer with the ability to control wall temperature. Note primary drying time.



Made controlled temperature walls. Also tried adding commercial insulation and aluminum shield to minimize edge effect.

Run ALPHA with a 19 vials and wall at +20C
Experimental Design
3ml fill in a 10 ml vial
Water
Freezing step ramp at 0.5C/min to -40C, hold 30 mins
Primary drying shelf set point 20C
Vacuum set point 100 mT

ALPHA with +20C Walls
532 minutes in primary drying. (reminder: target is 636 mins)

Run ALPHA with a 19 vials and wall at -20C
Experimental Design
3ml fill in a 10 ml vial
Water
Freezing step ramp at 0.5c/min to -40C, hold 30 mins
Primary drying shelf set point 20C
Vacuum set point 100 mT

ALPHA with -20C Walls
557 minutes in primary drying (reminder: target is 636 mins)

Experimental Group #2

5% Sucrose

Eliminate radiant energy to minimize edge effect vials by cooling and/or shielding the chamber walls of the freeze dryer.



Run ALPHA with a 19 vials and wall at +20C

Experimental Design
2ml fill in a 10 ml vial
5% sucrose by weight
Freezing step -10 hold 30 mins, -40C hold 120 mins
Primary drying shelf set point -20C
Vacuum set point 60 mT



ALPHA, wall at +20C
789 minutes in primary drying

Array and thermocouple placement in Alpha



Experimental Group #2

5% Sucrose Continued


Run ALPHA with a 19 vials and wall at -20C

Experimental Design
2ml fill in a 10 ml vial
5% sucrose by weight
Freezing step -10 hold 30 mins, -40C hold 120 mins
Primary drying shelf set point -20C
Vacuum set point 60 mT



ALPHA, wall at -20C
1096 minutes in primary drying

Results: Cooling the chamber wall to -20C appeared to make some difference but the elimination of radiant energy does not appear to be the answer to simulating a larger freeze dryer. Why?



Experimental Group #3

Water

Run systems with wall at -40C. Study the rate of sublimation in individual vials in a freeze dryer that has eliminated radiant energy to gain a better perspective on what is happening in the dryer.


Results

Low temperature chamber walls does not eliminate "edge effect"


Theory: Eliminating radiant energy did not appear to eliminate edge effect in the 19 vial array. Taking a closer look at how the individual vials are responding may help us understand this further. The common perception is that when we eliminate radiant energy the edge effect should be eliminated.

Third Experimental Set: Utilizing the gravimetric measurement method determine the rate of freeze drying by various vials within the array. Do perimeter vials freeze dry differently than center vials even minimal radiant energy is present?

Run ALPHA with a 19 vials and wall at -40C
Experimental Design
3ml fill in a 10 ml vial
Water
Freezing step at 0.5 to -40C, hold 30 mins
Primary drying shelf set point +20C
Vacuum set point 100 mT



Run ALPHA with a 19 vials and wall at -40C with aluminum shield and additional insulation. (10+ runs conducted)
Experimental Design
3ml fill in a 10 ml vial
Water
Freezing step at 0.5 to -40C, hold 30 mins
Primary drying shelf set point +20C
Vacuum set point 100 mT



% dry by weight 25% of the way into primary drying.

Results: Edge effect NOT eliminated by lowering the temperature of the chamber wall. If the presumption that radiant energy caused edge effect and radiant energy is eliminated, then edge effect should be eliminated. This experiment does not verify this assumption. Gravimetric studies indicate that edge vials are drying faster than the other vials, even when there is minimal radiant energy

Experimental Group #3

5% Sucrose

Run systems with wall at -40C. Study the rate of sublimation in individual vials in a freeze dryer that has eliminated radiant energy to gain a better perspective on what is happening in the dryer.

Results

Low temperature chamber walls does not eliminate "edge effect" as evidence by run time not approaching REVO run time of 2026 mins


Run ALPHA with a 19 vials and wall at -40C

Experimental Design
2ml fill in a 10 ml vial
5% sucrose by weight
Freezing step -10 hold 30 mins, -40C hold 120 mins

Primary drying shelf set point -20C
Vacuum set point 60 mT

Gravimetric tests not yet conducted.

ALPHA, wall at -40C
1042 minutes in primary drying




Experimental Group # 4

Water

Adjacent vials affect the thermal characteristics of one another. Since edge vials do not have the same number of adjacent vials they sublime differently. Have the controlled temperature wall come in contact with the vials to eliminate edge effect.

Theory: The product vials have a thermal effect on one another (essentially acting as cold sinks and hot sinks). Simulate adjacent vials by directly contacting the edge vials and controlling the temperature.


Fourth Experimental Set: Find a way to contact the cold temperature walls to the outside edge vials. If "edge effect" is eliminated this indicates that vial to vial contact has a significant effect in "edge effect".



No longer using the Alpha unit. Designed and build the MicroFD with "thermal simulator" temperature controlled blocks.

Run MicroFD with a 19 vials and wall at -40C

Experimental Design
3ml fill in a 10 ml vial
Water
Freezing step ramp 0.5 to -40C, hold 30 mins
Primary drying shelf set point +20C
Vacuum set point 100 mT



Std Dev 4.4

MicroFD with simulator blocks (aka LyoSim™)
633 minutes in primary drying
ONLY 3 MINUTES DIFFERENT THAN THE FULL TRAY RUN DONE IN THE REVO

EDGE EFFECT ELIMINATED
ALL VIALS NOW BEHAVE LIKE CENTER VIALS
The use of thermal simulator blocks in the MicroFD created an environment where the sublimation rate was very similar (essentially the same) in every vial. The elimination of edge effect indicates the MicroFD can be utilized for protocol development for larger freeze dryers based on its ability to simulate center vials.

Summary

Water

Experiment #	Product	Conditions (see experiments for details)	Time	Vial variation in sublimation rate
1	Water	REVO full tray	636 mins	Not studied
1	Water	REVO 19 vials	512 mins	High disparity, data not reported here
2	Water	Alpha 19 vials +20C walls	532 min	High disparity, data not reported here
2	Water	Alpha 19 vials -20C walls	557 mins	High disparity, data not reported here
3	Water	Alpha 19 with -40C wall	Not timed	High disparity, edge effect seen (see pg 3)
3	Water	Alpha 19 vials with -40C wall and additional radiant energy shielding	Not timed	High disparity, edge effect seen (see pg 3)
4	Water	MicroFD with LyoSim (thermal sink emulators touching vials)	633 mins (successfully simulates REVO run with full tray)	Variation in vials minimized, edge effect eliminated (see pg 4)

(see specific experiments for additional details)

Patent # 9121637 Accuflux Heat flux sensor in a freeze dryer
LyoSim Patent Pending

Summary

5% Sucrose

<i>Experiment #</i>	<i>Product</i>	<i>Conditions (see experiments for details)</i>	<i>Time</i>	<i>Vial variation in sublimation rate</i>
1	5% Sucrose	REVO full tray	2026 mins	Not yet studied
2	5% Sucrose	Alpha 19 vials +20C walls	789 mins	Not yet studied
2	5% Sucrose	Alpha 19 vials -20C walls	1096 mins	Not yet studied
3	5% Sucrose	Alpha 19 with -40C wall	1042 mins	Not yet studied
4	5% Sucrose	MicroFD with LyoSim (thermal sink emulators touching vials)	Not yet studied	Not yet studied

Product Experiments with customer

Additional & Subsequent Studies.

Big Pharma Collaboration: As part of our on-going work we conducted a blind study with a large pharma group. They created a cycle in a SP Scientific LyoStar. They supplied us with vials and product that had a critical temperature of -24C. Within 3 runs we were able to fine tune the Micro FD to simulate their LyoStar freeze dryer run to their satisfaction. They have subsequently purchased a unit for further study.