

The logo features the text "SiO2" in a large, white, sans-serif font. The "i" is lowercase, while "S" and "O" are uppercase. The "O" is a simple circle. The "2" is a simple numeral. Below this, the words "MATERIALS SCIENCE" are written in a smaller, blue, sans-serif font, followed by a small "SM" trademark symbol.

SiO₂

MATERIALS SCIENCESM

Next Generation Vials for COVID-19 vaccines and cell and gene therapies.

Millrock Webinar

Chris Weikart & Peter Sagona

April 8th, 2021

Challenges with COVID-19 Vaccines and Gene Therapy Drugs

- **Biologic Based Drugs/Vaccines**
 - Complex molecules
 - Inherently unstable in liquid formulations
- **Currently requires cryogenic storage to meet shelf life**
- **Efforts ongoing to Lyophilize (Freeze Drying) drug to eliminate cold storage.**

Technology

Hybrid Parenteral Vials – Materials of Construction



□ Primary Container (polymer)

- Injection Stretch Blow Molded
- Medical Grade Cyclic Olefin Polymer (COP/COC)
- Low dimensional Variability
- Optically clear & shatter resistant
- High barrier to water vapor

□ Barrier Coating System (glass-like)


- Plasma Enhanced Chemical Vapor Deposition
- Silicon-based coating
- High barrier to oxygen
- Inert No extractable/no leachables
- Barrier to label adhesives/polymer additives
- Optically clear




 Barrier Coating System
 Primary Container

Proprietary Technology Applied to Any Container Geometry

- Plasma enhanced chemical vapor deposition (PECVD)
- Coating applied to the inside of syringes, vials, cartridges of all geometries
- Inert drug contact surface is equivalent irrespective of container



1mL staked needle syringe



6mL vial

Characteristics of Barrier Coating System

- ❑ Chemically Inert & hydrophobic
 - Hydrolytically resistant (low Si dissolution)
 - pH resistant (no ion migration)
- ❑ Covalently bound to polymer
- ❑ Barrier to gases
 - Oxygen
 - Ethylene oxide




■ Barrier Coating System
■ Primary Container


“Covalent” Interfacial Bonds

- | | |
|--|---|
| <ul style="list-style-type: none"> • Si-O-Si • Si-CH_x | <ul style="list-style-type: none"> • O₃-Si-C • Si-O-C • Si-C • Si-CH_x |
|--|---|

“Hydrophobic” Surface Wettability



“Organosiloxane” Coating Chemistry



Lyophilized Drugs/Vaccines


Exploiting SiO₂ Hybrid Vials for Lyo

	Borosilicate Glass	SiO ₂ Medical Products
No Breakage	<ul style="list-style-type: none"> • Freezing expansion/contraction • Filling, handling & transport 	<ul style="list-style-type: none"> • Shatter resistant hybrid material • Breakage eliminated at high fill
Lower Particles	<ul style="list-style-type: none"> • Wall Shear from freeze expansion • Glass particles injected into patient 	<ul style="list-style-type: none"> • Covalent bonding • Lower shear force • Eliminates glass particles
Drying Consistency	<ul style="list-style-type: none"> • Variability in mass and dimensions • Higher heat transfer variability 	<ul style="list-style-type: none"> • Lower mass and dimensional variation • Flatter bottom vials
No Wall Residue	<ul style="list-style-type: none"> • Hydrophilic surface • Solution wicks & cake sticks to wall 	<ul style="list-style-type: none"> • Consistent hydrophobic surface • No cake sticking or wicking
Less Protein Aggregation	<ul style="list-style-type: none"> • Higher shear stress on proteins • Protein denaturing & aggregates 	<ul style="list-style-type: none"> • Lower protein shear forces • Lower aggregates & denaturing
Higher Fill Volumes	<ul style="list-style-type: none"> • Fill volumes kept below 50% • Risk of breakage 	<ul style="list-style-type: none"> • 100% fill volumes possible • No breakage

Polymer Thermal Conductivity and Vial Base Reduces Heat Transfer

Heat transfer coefficient of vials (K_v)
is dependent on:

- ❑ Thickness & mass of vial wall
- ❑ Thermal conductivity (h) of material
 - COP/COC – 0.13 W/m K
 - Glass – 1.20 W/m K
 - Air – 0.023 W/m K
- ❑ Contour of vial base



$$K_v \approx h/\text{thickness}$$


Improving COP/COC Vial Heat Transfer

Base of SiO₂ Vials

Standard



Flat-Bottom




Ink-blot
test for
flatness

Lower Dimensional and Mass Variability

Mass Variability


Vial Type	Mass (g)
Glass	11.708 ± 0.085
COP	6.726 ± 0.006
Coated COP	6.728 ± 0.005

Dimensional Variability



Improved Heat Transfer and Less Batch Variability


Tray Capacity (240 10ml Vials)



● Tested vials

Vial Type	$K_v \cdot 10^4$ (cal/s/cm ² /°C)	Standard Deviation
Glass	4.23	± 0.19
Coated COP	3.18	± 0.07
Coated COP (Flat Bottom)	3.56	± 0.07

Lower Moisture Loss Variability Within Batch



Relative Std Dev (% water lost)		
Vial Type	▼ 18 hrs	▼ 23 hrs
Glass	6.5	0.2
Coated COP (Flat Bottom)	1.7	0.5

- ❑ Glass has more water loss variability during primary drying (sublimation).
- ❑ Ensures more consistent lyo product quality.

No Breakage Irrespective of Fill Volume


Fill Volume (mL)	Fill Percentage (%)	Glass		SiO2	
		# Broken	%	# Broken	%
10	100	4/12	33	0/16	0
8	80	5/27	19	0/37	0
6	60	1/38	3	0/38	0
4	40	0/35	0	0/38	0



- ❑ Vials filled with 10% w/v mannitol solution
- ❑ Freezing: Shelves were cooled to -40°C at a rate of 1°C/min and held for two hours to freeze
- ❑ Drying: Chamber pressure was reduced to 150mT and the shelf was warmed to 5°C at 0.5°C/min.

NOTE: Conducted at Millrock Technology

Fewer Subvisible Particles after Reconstitution



*Particle analysis by FlowCAM® microflow digital imaging which detects particles between 2 -100 μm

Similarities and Improvements Compared to Glass


Lyo Formulation:

1 mg/ml IVIg
10mM glycine
5% w/v sucrose
0.02% v/v polysorbate 20



SiO₂ Coated

- ❑ **Similarities:** Initial residual moisture content, cake appearance, reconstitution time, and protein recovery.
- ❑ **Improvements:**
 - 4X Lower particle levels in reconstituted formulations in coated COP.
 - Less variable Kv and drying rates.



Eliminates Wall Residue – Hydrophobic Coating



SiO₂ Coated

Glass



❑ Placebo Formulation:

- hemophilia
- high salt content
- polar solution

❑ Wall Residue

- Silica-coated vials: NONE
- Glass: Significant

❑ Conclusion: hydrophobic coating surface reduces wall residue for some formulations.

Cold Storage

Exploiting SiO₂ Hybrid Vials for Cell & Gene Therapies

- ☐ No CCl Leakage at Cryo Temps (-196°C)
- ☐ Low Particles
- ☐ No Organic or Inorganic Leachables
- ☐ Barrier to Oxygen
- ☐ No Breakage
- ☐ Serialization – Track & Trace
- ☐ Ready-to-use


Freeze-Thaw Study Design (-20°C & -70°C)

Experimental:

- Vials filled with various volumes of high purity water and enclosed with rubber stopper and aluminum crimp
- Vials then placed in freezer, one group at -20°C and another group at -70°C
- After a 24-hour at temperature, vials were removed from freezer and held at room temp for 24 hours.
- Freeze-thaw cycle repeated 3 times with visual inspections after each cycle (cycle = low temp + room temp)




10 mL Vial Fill
(representative)



No Breakage Observed for SiO₂ Vials at -70°C

Low temp. to room temp. (24 hr. soak at each temp, no ramp between)

Manufacturer	Material	Fill Volume	# Failures (n = 50)			
			Cycle 1	Cycle 2	Cycle 3	TOTAL
-70C Cycling						
SiO2	COP 690R	6.5 mL	0	0	0	0
SiO2	COC S04	6.5 mL	0	0	0	0
Schott	Glass	6.5 mL	0	3	1	4
OMPI*	Glass	6.5 mL	2	1	2	5




No Breakage for any vials at -20°C cycling

*Ompi vials are 8 mL vial geometry

Deep Cold Storage CCI Experimental

CO₂ Headspace Analysis @ -80C




CCI Study Design (-80°C)

Set	Type	Conatiner Closure System (CCS)			Vial-Closure	Seal	Storage	Time Points							Total Meas
		Vial	Cap	Stopper	Finish	Compression	Temp.	t=0	1	2	1	3	6	12	
					(mm)	(RSF-lbf)	(C)		week	weeks	month	months	months	months	
1	Test Samples	SiOPlas 10 mL COP Vial	West Flip-Off Seal (Aluminum Crimp)	West NovaPure 1343, 4023/50	20	Low Mid High	-80	360	20	20	20	20	20	20	720
									20	20	20	20	20	20	
									20	20	20	20	20	20	
2	Test Samples	SiOPlas 10 mL COP Vial	West Flip-Off Seal (Aluminum Crimp)	West NovaPure S10-F451-4432/50	20	Low Mid High	-80	360	20	20	20	20	20	20	720
									20	20	20	20	20	20	
									20	20	20	20	20	20	
3	Test Samples	SiOPlas 10 mL COP Vial	West Flip-Off Seal (Aluminum Crimp)	West Westar RU (non-Flurotec coated) S-1727, 4432/50, B2-42	20	Low Mid High	-80	360	20	20	20	20	20	20	720
									20	20	20	20	20	20	
									20	20	20	20	20	20	
4	Positive Controls	SiOPlas 10 mL COP Vial (5µm Laser Drilled)	West Flip-Off Seal (Aluminum Crimp)	West NovaPure 1343, 4023/50	20	High	-80	30	5	5	5	5	5	5	60

- All test samples stoppered and crimped at Genesis Packaging Technologies under ambient conditions
- Positive controls samples, 5 µm laser drilled SiOPlas 10 mL COP vials, prepared by Lighthouse Instruments
- All headspace CO₂ measurements via FMS-Carbon Dioxide Headspace Analyzer performed by Lighthouse Instruments




No measurable changes in headspace observed after 6 months at **-80°C**




Cryogenic Storage CCI Experimental

O₂ Headspace Analysis @ **-180C**



Results: 10 mL Vial Headspace O₂ Concentration (%) after 3 months at **-180°C**

West NovaPure 4432 Stopper + Aluminum Crimp



Test Samples: n=15
Positive Controls: n=5



Conclusion: No measurable changes in headspace observed after 3 months at -180°C

Cell & Gene Therapy Company Using SiO₂ (-80°C Cold Storage)



- ❑ Peptide-based viral vector formulation
 - soluble adaptor to redirect human T cell engineered to express an universal chimeric antigen receptor UniCAR against PSMA expressing solid tumors
 - Packaged in 2mL SiO₂ Vials
 - 1+ year into 3-year cold storage (i.e. -80°C) stability study.
- ❑ Eliminated borosilicate glass as a packaging option due to problems with ion migration into drug product.

Demonstrated Drug Stability in SiO₂ 2mL Vial (1+ years @ -80°C Storage)



Purpose/Parameter	Method	Acceptance criteria	t0	1 m	3 m	6 m	9 m	12 m
See				□	□	□	□	□
Appearance	Ph. Eur. 2.1.1 / Ph. Eur. 2.2.2.	colorless and clear	colorless and clear	colorless and clear	colorless and clear	colorless and clear	colorless and clear	colorless and clear
Identity	¹ H-NMR	For information only	conforms	conforms	conforms	conforms	conforms	conforms
Purity	RP-HPLC	For information only	97.5*	97.4	97.5	97.1	97.6*	98.5*
Potency/identity	Standardized cellular flow assay	50-150 %	94	105	79	101	101	125
Content	UV (280 nm)	0.6 mg/mL ± 10 %	0.60	0.60	0.60	0.61	0.60	0.61
pH	Ph. Eur. 2.2.3.	4.5 ± 0.2	4.4	4.5	4.5	4.4	4.5	4.4
Sterility	Ph. Eur. 2.6.1.	Sterile, no microbial contamination	Sterile	-	-	-	-	-
Particulate contaminants (sub-visible particles)	Ph. Eur. 2.9.19.	≥ 10 µm: < 6000 per container ≥ 25 µm: < 600 per container	64 2	-	-	-	-	43/container 0/container
Endotoxins	Ph. Eur. 2.6.14	< 5 EU/mL	0.15 EU/mL	-	-	-	-	-
Test for integrity	According to ASTM E 515 – 2011	no leakage detected		-	-	-	-	-

Summary – SiO₂ Hybrid Vials

- ❑ Developed for Lyophilized Drug/Vaccine products
 - Lower heat transfer variability → Lower freezing/drying variability
 - No breakage and no wall residue
 - Fewer particles and aggregates compared to glass
- ❑ Demonstrated to withstand extreme cold storage requirements.
 - No CCI leakage down to -180°C
 - Drug stability with a viral vector formulation stored at -80C for 1+ years