

Water Vapor Pressure Table

Temp C	Press mT
0	4579
- 4	3280
- 8	2326
- 12	1632
- 16	1132
- 20	930
- 24	526
- 28	351
- 32	231
- 36	150
- 40	96.6
- 44	60.9
- 48	37.8
- 52	23.0
- 56	13.8
- 60	8.0
- 64	4.6
- 68	2.6
- 72	1.4
- 76	0.77
- 80	0.40
- 84	0.20
- 88	0.10

Courtesy of Millrock Technology, Inc